

# XCamShift

XCamShift (XCS) is a new implementation of the camshift chemical shift force field<sup>1</sup> for XPLOR-NIH<sup>2</sup>. The camshift implementation has a number of interesting and useful features

- It is very thoroughly tested against the original camshift implementation in almost<sup>3</sup> and has a large and full test suite which checks that XCS produces results that are identical with camshift in almost.
- XCS uses cython<sup>4</sup> to give a forcefield with native performance written in python
- XCS provides extensions against the standard camshift forcefield including the ability to carry out ensemble calculations, weighting of shift terms and variable well widths to allow for variation in statistical errors in measured chemical shift.
- XCS is modular and individual chemical shift terms can be re weighted during calculations and there is scope for whole components of the force field [e.g hydrogen bonding, aromatic ring currents etc] to be modified or replaced independently.
- XCS provides simple human readable files defining the components of the force field.
- XCS is open source and is free software under the Library Gnu Public License (LGPL)

For information on how to install, compile and use XCS see the sections below.

## The Details

### Installation

XCamShift can be installed from pre compiled binaries or can be compiled on the users own computer targeting a particular installation of xplor-nih. The

---

<sup>1</sup>**camshift:** Kohlhoff, K.J. et al. 2009. Fast and accurate predictions of protein NMR chemical shifts from interatomic distances. *Journal of the American Chemical Society*. **131** (39),pp.13894–13895. Robustelli, P. et al. 2010. Using NMR chemical shifts as structural restraints in molecular dynamics simulations of proteins. *Structure*. **18** (8),pp.923–933. [Kai Kohlhoff's thesis](#) "Protein Chemical Shifts as Structural Restraints in Molecular Dynamics Simulations, University of Cambridge, England, May, 2008

<sup>2</sup>**xplor-nih:** Schwieters, C.D. et al. 2003. The Xplor-NIH NMR molecular structure determination package. *J. Magn. Reson.* **160** (1),pp.65–73. Schwieters, C.D. et al. 2006. Using Xplor NIH for NMR molecular structure determination. *Progress in Nuclear Magnetic Resonance Spectroscopy*. **48** ,pp.47–62.

<sup>3</sup>**almost:** Fu, B. et al. 2014. ALMOST: An all atom molecular simulation toolkit for protein structure determination. *Journal of computational chemistry*. **35** (14),pp.1101–1105.

<sup>4</sup>**cython** C extensions for python (an optimising static compiler for both the Python programming language and the extended Cython programming language)

precompiled binary python modules includes all the software required to run XCS and you should not require any other software to be installed other than XPLOIR-NIH

1. OSX x86\_64 (64 bit)
2. OSX i686 (32 bit) - coming
3. Linux x86\_64 (64 bit)
4. Linux x86\_64 (32 bit) - coming
5. sources: see downloading the source from Github below

before you install you should

1. if md5sum is installed checks the md5sum of the download is correct md5sum XCamShift\_XXX.tgz [on OS X you will need md5 -r XCamShift\_XXX.tgz instead]. The md5sum is stored in XCamShift\_XXX\_md5.txt
2. extracts the archive XCamShift\_XXX.tgz into a convenient temporary directory using the command tar -zxvf XCamShift\_XXX.tgz

the easiest way to install XCamShift is to run the installation command

```
install_dist.sh <dist_dir> <install_dir>
```

which is included with the distribution. When using this command replace <dist\_dir> with the path to the directory that contains the xcamshift distribution (this is the directory that contains the install\_dist.sh file itself) and install\_dir is the directory that contains the particular xplor distribution to install into.

**NOTE 1.** if you want to use XCamShift without modifying your xplor-nih installation see the section [setting paths](#) below

---

**NOTE 2.** In summary this is what install\_dist.sh does ( in the following xxx refers to a specific version e.g. v1, <platform> is a computer platform e.g. OSX or Linux and <architecture> is a cpu architecture e.g. x86\_64 or i686)

1. move the .so files in xcamshift\_XXX/modules into the folder python/bin.<platform>\_<architecture> in you xplor installation; on OSX <platform> will be Darwin\_<XX> (where <XX> is the darwin version) and just Linux for linux. <architecture> will be either X86\_64 or i686 for 64 and 32 bit computers respectively.

2. Move the directory database/XCamShift which contains the forcefield definition into the database directory of your XPLORE-NIH installation.
3. copy the contents of the camshift\_XXX/python directory into the python directory of your xplor distribution
4. test the installation by running the script src/suite.py with the command line pyXplor src/suite.py. You shouldn't see any errors of the form fail or error or abnormal program termination; any of these show a bug is present and should be reported as described in bug reporting below.

## Using XCamShift

XCS can only be accessed from the python interface of xplor NIH currently. Here are some snippets showing how to carry out the main operations required to run XCS.

importing shifts

```
csFile=open("test_data/alvin/cs.dat")
cs=csFile.read()
csFile.close()
```

the format of a chemical shift restraint file is detail below under restraints  
setting up a calculation

```
from XCamShift import XCamShift
XCamShift.addRestrains(cs)
potList.append(XCamShift)
```

setting up weights

```
XCamShift.setScale(6) # this will change the weights of all terms
```

Format of the xcamshift restraint list

camshift restraint files are typical xplor restraints files and accept acclamation marks (!) as comments that run to the end of lines. The format is

```
<weight-statement>|<class-statement>| <assign-statement>
```

```
weight-statement = WEIGHt <float>
class-statement = CLASs <four-letter-code>
assign-statement = ASSIgn <selection> <float>|
                  ASSIgn <selection> <float> <float>
```

the weight and <class> statement set weights for all following <assign> statements. statements can take either one and two float values which are in order: the chemical shift and the error of the chemical shift. It is not possible to set individual weights without also setting an error without recourse to the weight statement. The class statement can be used to tag shifts for later manipulation (currently you can just use this to reset the weight of a set of restraints).

An example of an XCamShift restraint file would be

```
class TEST                                ! this is a comment
assign ( resid 2 and name C ) 177.477      ! restraint with default error
weight 2.0
assign ( resid 3 and name C ) 175.002 0.1    ! error is 0.1 weight 2.0
weight 0.5
assign ( resid 4 and name C ) 176.647 1.0    ! error 1.0 weight 0.5
```

It should be noted that chemical shift restraint files can't currently take sophisticated expressions as provided by the classic xplor syntax.

notes:

1. restraints from the first and last residues in a segment are ignored and a warning issued as these shift cannot be used by camshift as it calculates shifts for the atoms in the central residue of a triple of residues
2. segments with less than 3 residues are ignored and a warning is issued again because camshift as it calculates shifts for the atoms in the central residue of a triple of residues

#### Restraint objects

Restraint objects can be retrieved from the XCamShift forcefield using the method getRestraints() on the XCamShift potential object. Restraint objects have the following attributes

- obs - the observed shift (can be set)
- calcd - the last calculated shift, if shifts have been calculated
- diff - obs - calcd
- err - the err on the observed shift this will be the 1/2 width of the square well where no forces will be applied (can be set)
- weight - the weight of the shift (can be set)
- comment - any end of line comment present when the shift was read in
- name - the name of the atom being restrained (segid resid restype atom-name)
- calcd2 - calcd\*\*2
- obs2 - obs\*\*2

## Restraint statistics

the XCamShift force field object provides the following methods to retrieve information about the number of restraint and violations

- rms - the rms of diff between obs and calcd (see below for how a threshold is set)
- violations - the number of violations
- numRestrains - the number of restraints

As the range of chemical shifts is not uniformly distributed the value of the threshold [set using the attribute threshold of the forcefield] for restraint violations has to be treated differently and is used as a multiplier of the error. Therefore a restraint is violated when the following equality tests true

```
|diff| > threshold * err + err
```

## Listing predicted shifts

The chemical shifts predicted by XCamShift can be printed using the methods print\_shifts of the potential. The shifts are printed in the following format

SEGID	RESID	RES	HA	CA	HN	N	C	CB
AGB3	2	GLN	4.5501	54.2330	8.7379	121.0727	175.1000	29.2632
AGB3	3	TYR	5.0569	56.5513	9.5938	123.5669	174.7284	40.6070
AGB3	4	LYS	4.8969	54.5661	9.4434	121.8343	174.0538	35.2424
AGB3	5	LEU	4.6559	52.7604	9.2490	124.7154	174.3615	43.3438
AGB3	6	VAL	4.2273	60.9981	8.9896	124.2064	174.2570	32.7248
AGB3	7	ILE	4.2590	59.6785	9.0951	126.6597	174.7850	37.6970
AGB3	8	ASN	4.6711	51.7622	9.3382	125.9824	173.3525	38.0188
AGB3	9	GLY	4.0863	45.1157	8.6680	110.2235	174.1825	.
AGB3	10	LYS	3.9673	58.6702	9.3420	122.0362	176.9272	31.7829
AGB3	11	THR	4.2896	63.8501	8.1219	113.1026	174.1789	69.9168
...								
AGB3	54	VAL	4.4473	59.3198	8.6110	121.8589	173.1164	32.8183
AGB3	55	THR	4.8757	60.3260	8.6678	120.6908	173.4372	71.6735

note that missing chemical shifts which are structurally possible e.g. GLY CB are replaced by a . and that segids are wrapped with a || pair to make spaces visible.

## Rebuilding XCamShift

You should only need to recompile XCamShift if you want to make changes to the code in XCamShift, re-compile it on a system which isn't supported or use

a different compiler. If you only need to recompile with a different compiler or different target version of xplor-nih you won't need cython, you can just recompile the C code produced by cython.

**Pre requisites** To rebuild the camshift tables (it is unlikely that you will need to do this!) if you are running on an older version of python you may need to build argparse<sup>5</sup> and add it to the PYTHONPATH

**Downloading the source from Github** source code can be downloaded from the github master branch using the following link: locsmith/xcamshift/archive/master.zip to clone the complete repository use the following command `git clone https://github.com/locsmith/xcamshift.git`

### Rebuilding the C code

1. Firstly make sure you have a compiler that is compatible with GCC. On linux GCC and ICC have been tested. On OSX I have only used GCC from the homebrew project not Clang.
2. Make sure you have installed the source for your distribution of XPLOR-NIH
3. Download the source code for the version of python you want to target and build it. To test which version of python explore is using run `<xplor-nih directory>/bin/pyXplor test.py` with test.py containing `import sys ; print sys.version`. This will give you the version of python, typically the first line of output should be of the form 2.7.5 (default, Mar 9 2014, 22:15:05) and in this case the version is 2.7.5.
4. edit build\_all.py and change the variables
5. run build\_all.py

**Setting Paths** XCamShift can be tested without modifying xplor-nih by adding the following to your PYTHONPATH

`<XCamShift-directory>/src/cython:<XCamShift-directory>/XCamShift/src`

where is the directory that contains the XCamShift distribution

If it is not installed you will also need to install PyYAML<sup>6</sup> and add its installation directory to the PYTHONPATH variable

---

<sup>5</sup> **argparse** part of python since version 2.7, for earlier versions of python

<sup>6</sup> **pyyaml**: an implementation of the yaml markup language for python, you may also wish to install libyaml to get a much faster runtime

**Testing XCamShift** If XCamShift is installed or the paths are correctly set XCamShift can be tested using the following command

```
pyXplor <XCamShift-directory>/src/test/test_XCamShift_print_shifts.py
```

if they are not installed already you will also have to add python paths for the following python packages if you want to run the tests

1. nanotime<sup>7</sup>
2. unittest2<sup>8</sup>

**Reporting bugs** Bugs happen! XCamShift is well test but new code, please report bugs (or make requests for enhancements or even submit enhancements!) on the XCamShift github site:

The author can also be contacted at g dot s dot thompson at leeds dot ac dot uk

## Caveats

1. XCamShift calculates chemical shifts for hydrogen bonds but doesn't calculate any forces in the same manner as the original version. The hydrogen bond shift calculator can be disabled using the following code `<xcamshift-potential>.remove_named_sub_potential('HBOND')` where `<xcamshift-potential>` is the xcamshift potential instance.

## Tested versions

- **Targeted version of Camshift:** XCamShift is currently tested against version 1.35.0 of camshift which is implemented in almost 1.04.
- **Supported versions of XPLOR-NIH:** XCamShift currently supports version of XPLOR-NIH from 2.31 to 2.35 and the precompiled binaries are compiled against XPLOR-NIH 2.35. Please note that for some versions of xplor-nih before 2.35 you may need to recompile from source as the interface for the ensemble code has changed (this definitely included xplor-nih 2.33).
- **Compiler environment:** XCamShift binaries are compiled with GCC version 4.8.2 on OSX and Cython 0.20.2 (Note Cython is only required if the python code is modified, XCamShift can be re-compiled from the C++ source without Cython being present).
- **Targetted platforms:** XCamShift can be used on OSX and Linux under either a 32Bit or 64Bit intel environment (however 32bit OSX is not tested).

---

<sup>7</sup>**nanotime:** a nano second resolution timing service for python

<sup>8</sup>**unittest2:** a back port of the unit testing code from python 3+

### **Citing XCamShift**

There is currently no paper for XCamShift (I hope there will be one soon!). Until there is please cite: **XCamShift** G.S.Thompson, Astbury Centre for Structural Molecular Biology, University of Leeds, UK [<http://github.com/locsmith/xcamshift>]

### **References**